
Bedienungsanleitung • Mode d'emploi • Instruction

OptiMIC


Aufbau der Zoom-Optik
Mise en place de l'optique zoom
Mounting of the zoom-optics

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Pos. 1

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
Pos. 2

Video: Mounting of the zoom-optic Z1 / Z2 (Detent + Mod. II)


- A) Unscrew the screw (5) of the fixing clamp on the vertical holder of the base (4).
- B) Unscrew and take off the screw (2) from the zoom-optic.
- C) Introduce the zoom-optic in the vertical holder of the base (4) to the stop ring (1A or 1B, depending on the execution).
- D) The stroke-mark (3) has to be in front, tighten the screw (5).
- E) Screw in the screw (2) and tighten it a little.

Align the zoom-optic:

- F) Adjust the zoom on the minimum magnification.
- G) Align a test-lead point (Pos. 1) on the horizontal line of the screen. Begin preferably with the left side.
- H) Move the test-lead point to the right side (Pos. 1).
 - Is the test-lead point still exactly on the horizontal line (Pos.1), the zoom-optic is aligned correctly.
If the test-lead point is **above or below** the horizontal line (Pos. 2), make the following steps:

Adjust the fault:

- a) Unscrew the screw (5) a little, move the zoom-optic to the left or right side for adjust the horizontal line, tighten the screw (5).
- b) Repeat the points "G" and "H" till the test-lead point moves exactly on the horizontal line (Pos. 1).

Microscope: Mounting of the zoom-optic Z1 / Z2 (Detent + Mod. II)

- A) Unscrew the screw (5) of the fixing clamp on the vertical holder of the base (4).
- B) Unscrew and take off the screw (2) from the zoom-optic.
- C) Introduce the zoom-optic in the vertical holder of the base (4) to the stop ring (1A or 1B, depending on the execution).
- D) The stroke-mark (3) has to be in front, tighten the screw (5).
- E) Screw in the screw (2) and tighten it a little.

Align the zoom-optic:

- F) Adjust the zoom on the minimum magnification.
- G) Align a test-lead point (Pos. 1) on the horizontal line of the crosswire. Begin preferably with the left side.
- H) Move the test-lead point to the right side (Pos. 1).
 - Is the test-lead point still exactly on the horizontal line (Pos.1), the zoom-optic is aligned correctly.
If the test-lead point is **above or below** the horizontal line (Pos. 2), make the following steps:

Adjust **the fault**:

- c) Unscrew the screw (5) a little, move the zoom-optic to the left or right side for adjust the horizontal line, tighten the screw (5).
- d) Repeat the points "G" and "H" till the test-lead point moves exactly on the horizontal line (Pos. 1).

**For an eventual recentering of the goniometer, see the instruction
"Microscope with goniometer"**

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